

# IDYS<sup>®</sup>-ALIF Factsheet

## ALIF SYSTEM


SIMPLE



SAFE

ROBUST

## ALIF TECHNOLOGY



Improves OR efficiency

Strong anterior stabilization

Provide excellent stability

### X-Rays



# FEATURES & BENEFITS


## Cage/Plate Assembly

provides excellent stability upon insertion




## Angular Variability

allows for strong anterior fixation



## Low Plate Profile




Idys™-ALIF TiVac Cage			Height Anterior			
Footprint D X W	Lordosis Angle	H10	H12	H14	H16	H18
24x26	06°	24522006-S	-	-	-	-
	10°	-	24522010-S	-	-	-
	14°	-	-	24522014-S	-	-
26x32	06°	24520006-S	24520206-S	-	-	-
	10°	24520010-S	24520210-S	24520410-S	24520610-S	-
	14°	-	24520214-S	24520414-S	24520614-S	-
	18°	-	-	24520418-S	24520618-S	-
30x35	06°	24521006-S	24521206-S	-	-	-
	10°	24521010-S	24521210-S	24521410-S	24521610-S	-
	14°	-	24521214-S	24521414-S	24521614-S	24521814-S
	18°	-	-	24521418-S	24521618-S	24521818-S
30x40	06°	-	24523206-S	-	-	-
	10°	-	24523210-S	24523410-S	24523610-S	-
	14°	-	-	24523414-S	24523614-S	24523814-S
	18°	-	-	-	24523618-S	24523818-S

## Available Lordosis



Height of Associated Cage					
	HA10	HA12	HA14	HA16	HA18
Plate	24430010-S	24430012-S	24430014-S	24430016-S	24430018-S



Screws	L20	24010520-S
	L25	24010525-S
	L30	24010530-S
	L35	24010535-S

∅ 5.0 mm  
Self-tapping & self-drilling

## Clariance, Inc.

4809 N. Ravenswood Ave  
Suite 119 Chicago, IL 60640  
T: + 1 (773) 868-7041  
contact@clariance-spine.us

www.clariance-spine.com