

About Clariance

Founded in 2007 by Alain Tornier
(Former president and CEO of Tornier),
Clariance is a spinal device company
committed to designing, manufacturing
and marketing innovative solutions for
the treatment of spinal disorders.

Driven by surgeons' expertise, the company
provides advanced surgical applications
focused on fusion devices and minimally
invasive spine surgery. Improving patient
outcomes is our fundamental and permanent
concern. Clariance, with corporate Headquarters
in France and the US, is rapidly expanding
internationally by building relationships with the
scientific and business communities.



Clariance, Inc.

4809 N. Ravenswood Ave
Suite 119
Chicago, IL 60640
T: +1 773-868-7041



contact@clariance-spine.us

Clariance Headquarters

18 Rue Robespierre
Beaurains, France 62217
T: +33 3 21 16 12 15



contact@clariance-spine.com



www.clariance-spine.com



CLARIANCE

Idys™-ALIF



© 2018 Clariance, Inc. All rights reserved

24TRIEN



Cage/Plate Assembly

Cage/Plate assembly
provides excellent stability upon insertion




Low Plate Profile




Angle Variability

Angular Variability
allows for strong anterior fixation



Implant Range

Cage Sizes




24x26 mm x H 10 - 12 - 14 mm

26x32 mm x H 10 - 12 - 14 - 16 mm

30x35 mm x H 12 - 14 - 16 - 18 mm

30x40 mm x H 12 - 14 - 16 - 18 mm



Available Lordosis



18°



14°



10°



06°

Plate Sizes



IDYS-ALIF Plate H10 mm

IDYS-ALIF Plate H12 mm

IDYS-ALIF Plate H14 mm

IDYS-ALIF Plate H16 mm

IDYS-ALIF Plate H18 mm

One-step locking mechanism

Screw Sizes



IDYS-ALIF Screw L20 mm

IDYS-ALIF Screw L25 mm

IDYS-ALIF Screw L30 mm

IDYS-ALIF Screw L35 mm

Ø 5.0 mm

Self-tapping & self-drilling